

Product Datasheet

GTF2B Antibody (orb1268371)

Description

GTF2B Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This GTF2B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 42-71 amino acids from the N-terminal region of human GTF2B.

Target

GTF2B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[Q00403](#)

NCBI


[Q00403](#)

Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis of GTF2B using rab...



Western blot analysis in mouse kidney ti...