



Product Datasheet

NFYB Antibody (orb1268385)

Description	NFYB Antibody	
Species/Host	Rabbit	 Western blot analysis in A549 cell line ...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IF, WB	
Immunogen	This NFYB antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 17-46 amino acids from the N-terminal region of human NFYB.	 Fluorescent confocal image of HeLa cell ...
Target	NFYB	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	23 kDa	
Uniprot ID	P25208	
NCBI	P25208	
Dilution Range	For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	