

Product Datasheet

Tssk6 Antibody (orb1268407)

Description

Tssk6 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This Mouse Tssk6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 204-232 amino acids from the C-terminal region of mouse Tssk6.

Target

Tssk6

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW



30 kDa

Uniprot ID
Q925K9
NCBI
Q925K9
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Western blot analysis of Tssk6 using rab...



Western blot analysis in NCI-H292 cell l...



Western blot analysis in mouse kidney ti...