

Product Datasheet

OTUD7B Antibody (orb1268425)

Description

OTUD7B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This OTU7B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 757-784 amino acids from the C-terminal region of human OTU7B.

Target

OTUD7B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


93 kDa

Uniprot ID
[Q6GQQ9](#)
NCBI
[Q6GQQ9](#)
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis
in MDA-MB453 cell ...



OTU7B Antibody
immunohistochemistry
anal...