

Product Datasheet

Scyl2 Antibody (orb1268511)

Description

Scyl2 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This Mouse Scyl2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 585-613 amino acids from the Central region of mouse Scyl2.

Target

Scyl2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

103 kDa

Uniprot ID

[Q8CFE4](#)

NCBI


[Q8CFE4](#)

Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western Blot
at 1:1000
dilution +
mouse ...



Western blot
analysis in
NCI-H460 cell
l...