

Product Datasheet

PNKP Antibody (orb1268550)

Description

PNKP Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This PNKP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 109-138 amino acids from the N-terminal region of human PNKP.

Target

PNKP

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


57 kDa

Uniprot ID
[Q96T60](#)
NCBI
[Q96T60](#)
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in Jurkat cell lin...



Western blot analysis of PNKP using rabb...