

Product Datasheet

Pkmyt1 Antibody (orb1268572)

Description

Pkmyt1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This Mouse Pkmyt1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 114-143 amino acids from the N-terminal region of mouse Pkmyt1.

Target

Pkmyt1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

54 kDa

Uniprot ID

[Q9ESG9](#)

NCBI


[Q9ESG9](#)

Dilution Range


For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Western blot analysis in HL-60 cell line...



Western blot analysis of Pkmyt1 using ra...



Western blot analysis in mouse heart tis...