

Product Datasheet

CDH4 Antibody (orb1268616)

Description

CDH4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications

Immunogen

This CDH4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 175-203 amino acids from the N-terminal region of human CDH4.

Target

CDH4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

100 kDa

Uniprot ID

P55283

NCBI


P55283

Dilution Range


For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Western Blot
at 1:2000
dilution +
HepG2 ...



Formalin-
fixed and
paraffin-
embedded
hum...