

Product Datasheet

LYPLAL1 Antibody (orb1268658)

Description

LYPLAL1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This LYPL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 88-117 amino acids from the Central region of human LYPL1.

Target

LYPLAL1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

26 kDa

Uniprot ID

[Q5VWZ2](#)

NCBI

[Q5VWZ2](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


130
95

55
36
28

Western blot analysis
in mouse Neuro-2a ...



LYPL1 Antibody
immunohistochemistry
anal...



Flow cytometric
analysis of Neuro-2a
cel...