


Product Datasheet

CHST14 Antibody (orb1268661)

Antibody Description Description		CHST14 Antibody	 Western blot analysis in mouse stomach t...
Species/Host	Rabbit		
Reactivity	Mouse		
Conjugation	Unconjugated		
Tested Applications	WB		
Immunogen	This CHST14 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 178-205 amino acids from the Central region of human CHST14.		
Target	CHST14		
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.		
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		
Note	For research use only		
Application notes	For WB starting dilution is: 1:1000		
Isotype	Rabbit Ig		
Clonality	Polyclonal		
MW	43 kDa		
Uniprot ID	Q8NCH0		
NCBI	Q8NCH0		
Dilution Range	For WB starting dilution is: 1:1000		
Expiration Date	12 months from date of receipt.		