

Product Datasheet

EDN2 Antibody (orb1268720)

Description

EDN2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This EDN2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 107-136 amino acids from the C-terminal region of human EDN2.

Target

EDN2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

20 kDa

Uniprot ID

[P20800](#)

NCBI

[P20800](#)

Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.


K562

72



Western blot analysis in K562 cell line

...



Western blot analysis of EDN2 using rabb...