

Product Datasheet

TNFRSF11B Antibody (orb1268799)

Description

TNFRSF11B Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications
Immunogen

This TNFRSF11B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 243-271 amino acids from the Central region of human TNFRSF11B.

Target

TNFRSF11B

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

46 kDa

Uniprot ID

O00300

NCBI


O00300

Dilution Range


For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50 For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in NCI-H460 cell l...



Confocal
immunofluorescent
analysis of T...



TNFRSF11B Antibody
immunohistochemistry

...