

Product Datasheet

IL13 Antibody (orb1268800)

Description

IL13 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications
Immunogen

This IL13 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 118-146 amino acids from the C-terminal region of human IL13.

Target

IL13

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


16 kDa

Uniprot ID
P35225
NCBI
P35225
Dilution Range


For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in NCI-H460 cell l...



Confocal
immunofluorescent
analysis of l...



IL13 Antibody
immunohistochemistry
analy...