

Product Datasheet

GRID1 Antibody (orb1268857)

Description

GRID1 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This GRID1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 846-875 amino acids from the C-terminal region of human GRID1.

Target

GRID1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


112 kDa

Uniprot ID
Q9ULK0
NCBI
Q9ULK0
Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

m.bladder
250


Western blot analysis
in mouse bladder t...



GRID1 Antibody
immunohistochemistry
anal...