


Product Datasheet

TRIM45 Antibody (orb1268860)

Antibody Description Description		TRIM45 Antibody	 Western blot analysis in mouse heart tis...	
Species/Host		Rabbit		
Reactivity		Mouse		
Conjugation		Unconjugated		
Tested Applications		WB		
Immunogen		This TRIM45 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 389-417 amino acids from the C-terminal region of human TRIM45.		
Target		TRIM45		
Preservatives		Supplied in PBS with 0.09% (W/V) sodium azide.		
Form/Appearance		Liquid		
Concentration		batch dependent		
Storage		Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		
Note		For research use only		
Application notes		For WB starting dilution is: 1:1000		
Isotype		Rabbit Ig		
Clonality		Polyclonal		
MW		64 kDa		
Uniprot ID		Q9H8W5		
NCBI		Q9H8W5		
Dilution Range		For WB starting dilution is: 1:1000		
Expiration Date		12 months from date of receipt.		