



Product Datasheet

TUBGCP5 Antibody (orb1268901)

Description	TUBGCP5 Antibody	
Species/Host	Rabbit	 Western blot analysis in mouse NIH-3T3 c...
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This TUBGCP5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 172-201 amino acids from the N-terminal region of human TUBGCP5.	
Target	TUBGCP5	 TUBGCP5 Antibody immunohistochemistry an...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	118 kDa	
Uniprot ID	Q96RT8	
NCBI	Q96RT8	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	