


Product Datasheet

EIF2B3 antibody (orb126895)

Description	Rabbit Polyclonal to EIF2B3.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NR50
Dilution Range	WB: 1:500-1:3000, ELISA: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from J...