

Product Datasheet

SDF2L1 Antibody (orb1268965)

Description

SDF2L1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This SDF2L1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 160-189 amino acids from the C-terminal region of human SDF2L1.

Target

SDF2L1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

24 kDa

Uniprot ID

[Q9HCN8](#)

NCBI

[Q9HCN8](#)


Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


SK-BR-3



Western blot analysis
in SK-BR-3 cell li...



Western blot analysis
in mouse Neuro-2a ...



SDF2L1 Antibody
immunohistochemistry
ana...