

Product Datasheet

MAGED4 Antibody (orb1268983)

Description

MAGED4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This MAGED4B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 328-357 amino acids from the Central region of human MAGED4B.

Target

MAGED4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

81 kDa

Uniprot ID
[Q96JG8](#)
NCBI
[Q96JG8](#)
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

250


130

95

72

55

Western blot analysis
in HepG2 cell line...



MAGED4B Antibody
immunohistochemistry
an...