

Product Datasheet

CD1C Antibody (orb1268995)

Description

CD1C Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This CD1C antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 226-253 amino acids from the C-terminal region of human CD1C.

Target

CD1C

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

38 kDa


Uniprot ID
P29017
NCBI
P29017
Dilution Range

For WB starting dilution is: 1:1000


Expiration Date

12 months from date of receipt.

NCI-H460



Western blot analysis in NCI-H460 cell l...



Western blot analysis of CD1C using rabb...