


Product Datasheet

Q91ZF4 Antibody (orb1268996)

Description	Q91ZF4 Antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	This cathepsin Q antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 192-221 amino acids from human cathepsin Q.
Target	Q91ZF4
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q91ZF4
NCBI	Q91ZF4
Dilution Range	For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis in mouse spleen ti...