

Product Datasheet

TRIM38 Antibody (orb1269001)

Description

TRIM38 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This TRIM38 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 423-451 amino acids from the C-terminal region of human TRIM38.

Target

TRIM38

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


53 kDa


Uniprot ID
000635
NCBI
000635
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.


 Western blot analysis
in HeLa cell line ...

 TRIM38 Antibody
immunohistochemistry
ana...

 Flow cytometric
analysis of HeLa cells
(...