



Product Datasheet

RPL3 Antibody (orb1269002)

Description	RPL3 Antibody	
Species/Host	Rabbit	 Western blot analysis in MDA-MB435 cell ...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This RPL3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 224-251 amino acids from the Central region of human RPL3.	
Target	RPL3	 RPL3 Antibody immunohistochemistry analy...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	46 kDa	
Uniprot ID	P39023	
NCBI	P39023	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	