



Product Datasheet

THNSL1 Antibody (orb1269040)

Description	THNSL1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in K562 cell line ...
Immunogen	This THNSL1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 694-722 amino acids from the C-terminal region of human THNSL1.	
Target	THNSL1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	83 kDa	
Uniprot ID	Q8IYQ7	
NCBI	Q8IYQ7	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	