

Product Datasheet

NETO2 Antibody (orb1269043)

Description

NETO2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This NETO2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 174-203 amino acids from the N-terminal region of human NETO2.

Target

NETO2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

59 kDa

Uniprot ID

[Q8NC67](#)

NCBI


[Q8NC67](#)

Dilution Range

For FACS starting dilution is: 1:25For WB starting dilution is: 1:2000For IHC-P starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.




Overlay histogram showing Jurkat cells

S...



Western Blot at 1:2000 dilution + Jurkat...



NETO2 Antibody (N-term) immunohistochemi...