

## Product Datasheet

### TASP1 Antibody (orb1269051)

## Description

TASP1 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested Applications

IHC-P, WB

### Immunogen

This TASP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 364-395 amino acids from the C-terminal region of human TASP1.

### Target

TASP1

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

44 kDa

### Uniprot ID

[Q9H6P5](#)

### NCBI


[Q9H6P5](#)

### Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

### Expiration Date

12 months from date of receipt.



Western blot analysis of TASP1 C-term po...



Formalin-fixed and paraffin-embedded hum...