



Product Datasheet

ERO1L Antibody (orb1269053)

Description	ERO1L Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This ERO1L antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 350-379 amino acids from the C-terminal region of human ERO1L.	Western blot analysis in CEM cell line l...
Target	ERO1L	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	ERO1L Antibody immunohistochemistry anal...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	54 kDa	
Uniprot ID	Q96HE7	
NCBI	Q96HE7	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	