

## Product Datasheet

### TRIM11 Antibody (orb1269062)

## Description

TRIM11 Antibody

## Species/Host

Rabbit

## Reactivity

Human

## Conjugation

Unconjugated

## Tested

IHC-P, WB

## Applications

## Immunogen

This TRIM11 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 300-329 amino acids from the C-terminal region of human TRIM11.

## Target

TRIM11

## Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

## Form/Appearance

Liquid

## Concentration

batch dependent

## Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

## Note

For research use only

## Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

## Isotype

Rabbit Ig

## Clonality

Polyclonal

## MW

53 kDa

## Uniprot ID

[Q96F44](#)

## NCBI


[Q96F44](#)

## Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

## Expiration Date

12 months from date of receipt.



Western blot analysis  
in HepG2 cell line...



TRIM11 Antibody  
immunohistochemistry  
ana...