

Product Datasheet

ARFIP1 Antibody (orb1269245)

Description

ARFIP1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This ARFIP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 147-176 amino acids from the Central region of human ARFIP1.

Target

ARFIP1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


42 kDa

Uniprot ID
P53367
NCBI
P53367
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis
in MCF-7 cell line...



ARFIP1 Antibody
immunohistochemistry
ana...