



## Product Datasheet

### RNF180 Antibody (orb1269321)

| Description         | RNF180 Antibody                                                                                                                                                                                                        |  |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Species/Host        | Rabbit                                                                                                                                                                                                                 | Western blot analysis<br>in Y79 cell line l...                                      |
| Reactivity          | Human                                                                                                                                                                                                                  |                                                                                     |
| Conjugation         | Unconjugated                                                                                                                                                                                                           |                                                                                     |
| Tested Applications | IHC-P, WB                                                                                                                                                                                                              |                                                                                     |
| Immunogen           | This RNF180 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 401-430 amino acids from the C-terminal region of human RNF180.                                               |  |
| Target              | RNF180                                                                                                                                                                                                                 |                                                                                     |
| Preservatives       | Supplied in PBS with 0.09% (W/V) sodium azide.                                                                                                                                                                         |                                                                                     |
| Form/Appearance     | Liquid                                                                                                                                                                                                                 |                                                                                     |
| Concentration       | batch dependent                                                                                                                                                                                                        | RNF180 Antibody<br>immunohistochemistry<br>ana...                                   |
| Storage             | Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures. |                                                                                     |
| Note                | For research use only                                                                                                                                                                                                  |                                                                                     |
| Application notes   | For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50                                                                                                                                             |                                                                                     |
| Isotype             | Rabbit Ig                                                                                                                                                                                                              |                                                                                     |
| Clonality           | Polyclonal                                                                                                                                                                                                             |                                                                                     |
| MW                  | 68 kDa                                                                                                                                                                                                                 |                                                                                     |
| Uniprot ID          | <a href="#">Q86T96</a>                                                                                                                                                                                                 |                                                                                     |
| NCBI                | <a href="#">Q86T96</a>                                                                                                                                                                                                 |                                                                                     |
| Dilution Range      | For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50                                                                                                                                             |                                                                                     |
| Expiration Date     | 12 months from date of receipt.                                                                                                                                                                                        |                                                                                     |