

Product Datasheet

DOC2A Antibody (orb1269344)

Description

DOC2A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This DOC2A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 277-306 amino acids from the Central region of human DOC2A.

Target

DOC2A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:10~50For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

44 kDa

Uniprot ID

[Q14183](#)

NCBI


[Q14183](#)

Dilution Range


For FACS starting dilution is: 1:10~50For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Flow cytometric analysis of K562 cells (...)



Western blot analysis in K562 cell line ...



DOC2A Antibody immunohistochemistry anal...