

Product Datasheet

COPS7A Antibody (orb1269367)

Description

COPS7A Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This COPS7A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 237-266 amino acids from the C-terminal region of human COPS7A.

Target

COPS7A

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


30 kDa

Uniprot ID
[Q9UBW8](#)
NCBI
[Q9UBW8](#)
Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

h.placenta
95


Western blot analysis
in human placenta ...



COPS7A Antibody
immunohistochemistry
ana...