

Product Datasheet

PAGE5 Antibody (orb1269379)

Description

PAGE5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This PAGE5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 102-130 amino acids from the C-terminal region of human PAGE5.

Target

PAGE5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

14 kDa

Uniprot ID

[Q96GU1](#)

NCBI

[Q96GU1](#)


Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.

NCI-H292



Western blot analysis
in NCI-H292 cell l...



PAGE5 Antibody
immunohistochemistry
anal...