

Product Datasheet

WHSC1L1 Antibody (orb1269401)

Description

WHSC1L1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This WHSC1L1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 188-217 amino acids from the N-terminal region of human WHSC1L1.

Target

WHSC1L1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

162 kDa

Uniprot ID

[Q9BZ95](#)

NCBI


[Q9BZ95](#)

Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western Blot
at 1:2000
dilution +
mouse ...



Western blot
analysis of
WHSC1L1
using r...