

Product Datasheet

SIX5 Antibody (orb1269411)

Description

SIX5 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This SIX5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 257-286 amino acids from the Central region of human SIX5.

Target

SIX5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

75 kDa

Uniprot ID

[Q8N196](#)

NCBI

[Q8N196](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For IF starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

95


55

36


28

17

Western blot analysis
in MDA-MB231 cell ...



SIX5 Antibody
immunohistochemistry
analy...



Fluorescent confocal
image of HeLa cell ...