

Product Datasheet

KRT12 Antibody (orb1269415)

Description

KRT12 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This KRT12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 442-471 amino acids from the C-terminal region of human KRT12.

Target

KRT12

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

54 kDa

Uniprot ID

[Q99456](#)

NCBI


[Q99456](#)

Dilution Range


For FACS starting dilution is: 1:25For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Overlay histogram showing K562 cells sta...



Western Blot at 1:2000 dilution + rat ey...



Western Blot at 1:2000 dilution + mouse ...