

Product Datasheet

GSTM5 Antibody (orb1269421)

Description
GSTM5 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This GSTM5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 14-43 amino acids from the N-terminal region of human GSTM5.

Target

GSTM5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


26 kDa

Uniprot ID
P46439
NCBI
P46439
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
of lysate from hum...



GSTM5 Antibody
immunohistochemistry
anal...



Flow cytometric
analysis of K562 cells
(...)