

Product Datasheet

CTTNBP2NL Antibody (orb1269423)

Description

CTTNBP2NL Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This CTTNBP2NL antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 135-164 amino acids from the N-terminal region of human CTTNBP2NL.

Target

CTTNBP2NL

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

70 kDa

Uniprot ID

Q9P2B4

NCBI


Q9P2B4

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot at
1:1000 dilution + ZR-
75...



CTTNBP2NL Antibody
immunohistochemistry



Western blot analysis
in MDA-MB453, ZR-7...