

## Product Datasheet

### ZP4 Antibody (orb1269424)

## Description

## ZP4 Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This ZP4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 445-474 amino acids from the C-terminal region of human ZP4.

### Target

ZP4

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

59 kDa

### Uniprot ID

[Q12836](#)

### NCBI


[Q12836](#)

### Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

### Expiration Date

12 months from date of receipt.



Western Blot at  
1:2000 dilution +  
NCCIT ...



Western blot analysis  
in HepG2 cell line...



ZP4 Antibdy  
immunohistochemistry  
analysi...