

Product Datasheet

FBLIM1 Antibody (orb1269428)

Description
FBLIM1 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Western blot analysis
in A549 cell line ...

Immunogen

This FBLIM1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 294-323 amino acids from the C-terminal region of human FBLIM1.

Target

FBLIM1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P
starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


41 kDa

Uniprot ID
Q8WUP2
NCBI
Q8WUP2
Dilution Range

For WB starting dilution is: 1:1000For IHC-P
starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



FBLIM1 Antibody
immunohistochemistry
ana...