

Product Datasheet

STAMBP Antibody (orb1269440)

Description

STAMBP Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This STAMBP antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 355-383 amino acids from the C-terminal region of human STAMBP.

Target

STAMBP

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


48 kDa

Uniprot ID
095630
NCBI
095630
Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis
in A549 cell line ...



Western blot analysis
in mouse cerebellu...



STAMBP Antibody
immunohistochemistry
ana...