


Product Datasheet

OR6N2 Antibody (orb1269442)

Description	OR6N2 Antibody	 h.placenta 130 12 55 36 28 Western blot analysis in human placenta ...
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This OR6N2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 244-273 amino acids from the C-terminal region of human OR6N2.	
Target	OR6N2	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	36 kDa	
Uniprot ID	Q8NGY6	
NCBI	Q8NGY6	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	