

Product Datasheet

MYOD1 Antibody (orb1269465)

Description

MYOD1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, WB

Applications

Immunogen

This MyoD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 211-240 amino acids from the Central region of human MyoD1.

Target

MYOD1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

[P15172](#)

NCBI


[P15172](#)

Dilution Range


For WB starting dilution is: 1:1000 For IF starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot
at 1:1000
dilution Lane
1: ...



Western blot
analysis in
Hela cell line
...



Fluorescent
confocal
image of
Hela cell ...