


Product Datasheet

COMT antibody (orb12695)

Description	Rabbit polyclonal antibody to COMT
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant fusion protein of human COMT(NP_009294.1).
Preservatives	PBS with 0.02% sodium azide,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	25kDa, 30kDa
Uniprot ID	P21964
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF/ICC: 1:50-1:200, IP: 1:20-1:50, FACS: 1:20-1:50
Expiration Date	12 months from date of receipt.



Western blot
analysis of
extracts of
var...