

Product Datasheet

PREX1 Antibody (orb1269514)

Description

PREX1 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This PREX1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1134-1163 amino acids from the C-terminal region of human PREX1.

Target

PREX1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

186 kDa

Uniprot ID

[Q8TCU6](#)

NCBI


[Q8TCU6](#)

Dilution Range


For FACS starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:1000

Expiration Date


12 months from date of receipt.



Overlay histogram showing U-2OS cells st...



Antibody staining PREX1 in human brain t...



Western Blot at 1:2000 dilution Lane 1: ...