

Product Datasheet

FCN3 Antibody (orb1269521)

Description

FCN3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This FCN3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 214-243 amino acids from the C-terminal region of human FCN3.

Target

FCN3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


33 kDa

Uniprot ID
075636
NCBI
075636
Dilution Range


For WB starting dilution is: 1:500For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western Blot at 1:500 dilution + Jurkat ...



Western blot analysis in K562 cell line ...



FCN3 Antibody immunohistochemistry analy...