

Product Datasheet

DAZ4 Antibody (orb1269526)

Description

DAZ4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This DAZ4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 90-119 amino acids from the Central region of human DAZ4.

Target

DAZ4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

65 kDa

Uniprot ID
Q86SG3
NCBI
Q86SG3
Dilution Range

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.

95

55

36

28

17

Western blot analysis
in A549 cell line ...



DAZ4 Antibody
immunohistochemistry
analy...