

Product Datasheet

KIAA1324 Antibody (orb1269529)

Description

KIAA1324 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This K1324 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 704-732 amino acids from the C-terminal region of human K1324.

Target

KIAA1324

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

111 kDa

Uniprot ID

[Q6UXG2](#)

NCBI


[Q6UXG2](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis
in MDA-MB435 cell ...



Western blot analysis
in mouse spleen ti...



K1324 Antibody
immunohistochemistry
anal...