



Product Datasheet

GAA Antibody (orb1269532)


Description

GAA Antibody


Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IHC-P, WB
Immunogen	This GAA antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 174-203 amino acids from the N-terminal region of human GAA.
Target	GAA
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For IHC starting dilution is: 1:25For WB starting dilution is: 1:1000
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	105 kDa
Uniprot ID	P10253
NCBI	P10253
Dilution Range	For IHC starting dilution is: 1:25For WB starting dilution is: 1:1000
Expiration Date	12 months from date of receipt.



Antibody staining GAA in Human liver tis...



Western blot analysis of lysate from Jur...



Western blot analysis of lysates from A5...