

Product Datasheet

TSSK6 Antibody (orb1269584)

Description

TSSK6 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications
Immunogen

This TSSK6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 198-228 amino acids from the C-terminal region of human TSSK6.

Target

TSSK6

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


30 kDa

Uniprot ID
Q9BXA6
NCBI
Q9BXA6
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of TSSK6 using rab...



TSSK6 Antibody immunohistochemistry anal...



Flow cytometric analysis of 293 cells (r...